


Stool sample collection instructions

Your doctor ordered a stool test. Please make sure you understand which test(s) your doctor ordered for you.

- Read all instructions before starting
- Keep the materials out of reach of children and pets
- Make sure the materials haven't expired and are in good condition. If needed, contact your doctor's office for a replacement
- If your sample collection container contains liquid:
 - Don't drink the liquid. If the liquid gets in your eyes, then wash with water and seek medical attention right away
 - Don't pour the liquid out—the liquid is needed for testing. Don't add any liquid
 - If the liquid spills, overfills, or gets contaminated, contact your doctor's office for a replacement
- If you're collecting samples for more than 1 test, use a separate vial or cup for each test

For best results, follow the steps below.

Sterile screw-top cup (color or design may be different)



Stool testing options

- Calprotectin (16796)
- *Clostridium difficile* (*C diff*) Toxin/GDH with Reflex to PCR (91664)
- *Helicobacter pylori* (*H pylori*) Antigen, EIA (34838)
- Lactoferrin, Qualitative (10156)
- Lactoferrin, Quantitative (17321)*
- Pancreatic Elastase-1 (14693)

Get ready

- Wash and dry your hands
- Label the cup (**name, date of birth, collection date and time**)

Collect

- Lift the toilet seat. Place the collection hat on the toilet bowl. Lower the toilet seat
- Pass stool into the collection hat. Keep urine and toilet water away from the stool

Transfer

- Twist the lid off the cup
- Use the provided wooden stick to transfer a walnut-sized amount of stool from the collection hat into the cup. For *C diff*, transfer 1 teaspoon of loose/watery stool into the cup
- Screw the lid onto the cup tightly

Finish up

- Place the cup in the zip-top bag and seal the bag tightly (don't remove the sheet from the bag)
- Empty and throw away the collection hat
- Wash and dry your hands

Store and return

- Refrigerate or freeze the sealed bag
- Return the sealed bag with your sample(s) to the nearest Quest Patient Service Center or your doctor's office **within 48 hours of collection**

Total-Fix[™] vial



Stool testing options

- *Cryptosporidium* Antigen, EIA (90389)
- *Cyclospora* and *Isospora* Examination (10018)
- *Giardia* and *Cryptosporidium* Antigen Panel (39441)
- *Giardia* Antigen, EIA (8625)
- Ova and Parasites, Concentrate and Permanent Smear (681)
- Ova and Parasites with *Giardia* Antigen, EIA (1748)

Get ready

- Wash and dry your hands
- Label the vial (**name, date of birth, collection date and time**)

Collect

- Lift the toilet seat. Place the collection hat on the toilet bowl. Lower the toilet seat
- Pass stool into the collection hat. Keep urine and toilet water away from the stool

Transfer

- Twist the lid off the vial
- Use the spoon under the lid to transfer small amounts of stool from the collection hat into the vial. Stop when the liquid in the vial reaches the **FILL LINE**
- Screw the lid onto the vial tightly


Finish up

- Place the vial in the zip-top bag and seal the bag tightly (don't remove the sheet from the bag)
- Empty and throw away the collection hat
- Wash and dry your hands

Store and return

- Keep the sealed bag at room temperature. Don't refrigerate or freeze
- Return the sealed bag with your sample(s) to the nearest Quest Patient Service Center or your doctor's office **within 48 hours of collection**

*Note: Lactoferrin, Quantitative (17321) is not currently available in New York.



Cary-Blair vial



Stool testing options

- Gastrointestinal Pathogen Panel, Real-Time PCR (38470)
- *Salmonella/Shigella* Culture, *Campy* EIA and Shiga Toxin Reflex to *E coli* O157 Culture (32114)

Get ready

- Wash and dry your hands
- Label the vial (**name, date of birth, collection date and time**)

Collect

- Lift the toilet seat. Place the collection hat on the toilet bowl. Lower the toilet seat
- Pass stool into the collection hat. Keep urine and toilet water away from the stool

Transfer

- Twist the lid off the vial
- Use the spoon under the lid to transfer small amounts of stool from the collection hat into the vial. Stop when the liquid in the vial reaches the **FILL LINE**
- Screw the lid onto the vial tightly

Finish up

- Place the vial in the zip-top bag and seal the bag tightly (don't remove the sheet from the bag)
- Empty and throw away the collection hat
- Wash and dry your hands

Store and return

- Keep the sealed bag at room temperature. Don't refrigerate or freeze
- Return the sealed bag with your sample(s) to the nearest Quest Patient Service Center or your doctor's office **within 48 hours of collection**



Important! Samples must be returned within 48 hours of collection.



Find a location, schedule appointments, track results, and build your health history—all in one place at **MyQuest.QuestDiagnostics.com** or the **MyQuest® app**.

QuestDiagnostics.com

Quest Diagnostics Incorporated and its subsidiaries (Quest) complies with applicable federal civil rights laws and does not discriminate on the basis of race, color, national origin, age, disability, or sex. ATTENTION: If you speak English, language assistance services, free of charge, are available to you. Call 1.844.698.1022. ATENCIÓN: Si habla español (Spanish), tiene a su disposición servicios gratuitos de asistencia lingüística. Llame al 1.844.698.1022. 注意: 如果您使用繁體中文 (Chinese): 您可以免費獲得語言援助服務。請致電 1.844.698.1022.

Quest® is the brand name used for services offered by Quest Diagnostics Incorporated and its affiliated companies. Quest Diagnostics Incorporated and certain affiliates are CLIA-certified laboratories that provide HIPAA-covered services. Other affiliates operated under the Quest® brand, such as Quest Consumer Inc., do not provide HIPAA-covered services.

Quest®, Quest Diagnostics®, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. All third-party marks—® and ™—are the property of their respective owners. © 2025 Quest Diagnostics Incorporated. All rights reserved. SC14265CON 12/2025